

# Shower Trolley Service Record

<b>HEALTH UNIT NAME</b>	
<b>PRODUCT SERIAL NUMBER</b>	
<b>ASSET NUMBER</b>	
<b>DATE OF LAST SERVICE</b>	
<b>SERVICE PERSON</b>	

SERVICE CHECK POINTS	PASS	FAIL	COMMENTS
<b>BATH LINER AND PILLOW</b>			
Is the liner torn or punctured			
Is the foam firm and intact			
Is the liner clean underneath			
<b>DRAIN HOSE</b>			
Is the hose clean and free of damage			
<b>CORNER BUMPERS</b>			
Are they worn or gouged			
<b>TROLLEY SIDE / END RAILS</b>			
Are there any end rails bent or broken			
<b>DROPSIDE OPERATION</b>			
Are any of the dropside pivot bushes worn			
Are all the catches locking in safely			
<b>CASTORS</b>			
Is any castor visibly damaged			
Do the trolley brakes work correctly			
Does the trolley steer wheel lock in			
<b>HIGH / LOW OPERATION</b>			
Is the trolley raising and lowering acceptable			
Are the hydraulics or electrics trouble free			
<b>BOLTS PINS AND SPRINGS</b>			
Are all bolts nuts and pins tight and in place			
Are all springs returning properly			
Have all moving parts been lubricated with grease			
<b>ACCESSORIES</b>			
Are any extra options damaged			

